


56663

Upright Dustpan and Broom Set, 330 mm, Medium, Blue



The contact lip is designed for easy and efficient sweeping into the dustpan. The outer front surface has a raised area to prevent debris from falling out of the dustpan. The enlarged rear and side walls enable the dustpan to hold larger volumes of debris. The very effective, lightweight broom can be used for sweeping both dry and wet debris. The brush pattern is designed using 2 filament diameters. The outer row of filaments is softer than the middle row, which gives a very effective sweeping quality. The filaments under the scraping edge are angled to facilitate sweeping close to walls and corners as well as under workbenches etc. An effective scraping edge has been implemented in the design to enable loosening stubborn debris from the surface.

Technical Data

Item Number	56663
Bristle stiffness	Medium
Material	Polypropylene Anodised Aluminium Polyamide Polyester
Complies with (EC) 1935/2004 on food contact materials¹	Yes
Complies with EU Regulation 2023/2006/EC of Good Manufacturing Practice	Yes
Complies with FDA Regulationl CFR 21¹	Yes
Complies with UK 2019 No. 704 on food contact materials	Yes
Complies with REACH Regulation (EC) No. 1907/2006	Yes
Complies with California Proposition 65	Yes
Complies with Halal and Kosher	Yes
PFAS, Phthalates and BPA intentionally added	No
Box Quantity	10 Pcs.
Quantity per Pallet (80 x 120 x 180-200 cm)	120 Pcs
Quantity Per Layer (Pallet)	20 Pcs.
Box Length	1000 mm
Box Width	410 mm
Box Height	260 mm
Length	330 mm
Width	200 mm
Height	88 mm
Net Weight	0.82 kg
Weight bag (Recycling Symbol "4" Low Density Polyethylene (LDPE))	0.06 kg
Weight cardboard (Recycling symbol "20" PAP)	0.0668 kg
Tare total	0.1268 kg
Gross Weight	0.95 kg
Cubic metre	0.005808 M3
Recommended sterilisation temperature (Autoclave)	121 °C
Max. cleaning temperature (Dishwasher)	93 °C
Max usage temperature (food contact)	80 °C
Max usage temperature (non food contact)	100 °C
Min. usage temperature³	-20 °C
Gtin-13 Number	5705022044524
GTIN-14 Number (Box quantity)	15705028044532
Customs Tariff No.	96039099

New equipment should be cleaned, disinfected, sterilised and any labels removed, as appropriate to its intended use, e.g. high risk vs. low risk food production areas, general hospital areas vs. intensive care units, before use.

1. See Declaration of Compliance for further details on food contact
3. Do not store the product below 0° Celsius.